NOVOX[®]

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RIMOZIONE DEL DISPOSITIVO
 Chose the NovoX® Andro size suitable for the penis size, using the size selection guide on the package. Make sure that the device covers the entire lesion area. For surgical wounds, NovoX® Andro should be applied immediately after suturing the lesion.
 Users should wash their hands thoroughly before and after use. Moreover, wear disposable gloves and use tools that are preferably sterile or otherwise subject to careful disinfection when applying the product.
 • Before applying the device, clean the lesion area as directed by health care personnel using physiological saline, sterile water or other appropriate solution. Pat dry with sterile gauze to remove excess water. Wipe off excess disinfectant if used.
 • Open the package containing the device by pulling the flaps at the markings on the package, immediately before using the device.
 • Take out the device, taking care not to touch the impregnated inside that will be in contact with the wounds.
 • The device can be cut to reduce its length if necessary, first thoroughly disinfecting the equipment used (or using sterile equipment) and taking care not to damage the substrate or the textile coverings on the inside and outside. The device's closed tip can also be cut off, if necessary in the judgement of the doctor. It is recommended to follow the doctor's instructions carefully before cutting the device. It is recommended not to cut the device unless explicitly instructed by your doctor.
 • Fit the device, taking care not to deform or crumple it excessively and ensuring that it sticks firmly to the area to be treated.
 • Wear underwear, taking care to use a panty liner to prevent any small spills from soiling the underwear.
 • As an alternative to the above, if catheterisation is present that makes it impossible to insert the device, the following steps can be taken after thorough disinfection of the tools used to cut the device (or using sterile tools):
 - Cut off the distal end to allow passage of the catheter.
 - Open the device longitudinally by cutting on one of the welds.
 - Apply the device to the penis by wrapping the device over its entire length.
 - Restore the tubular shape of the device with at least 3 or more stitches (depending on the length) so that the device is well adhered to the penis.
 - For the execution of the suture referred to in the previous point, it is recommended to:
 - Use detached stitches (minimum 3)
 - Use slowly absorbable (e.g. polyglyconate) or non-absorbable thread of the same size used for tissue suturing.
DEVICE REMOVAL
 • To remove the device, simply pull it off. If removal is difficult, saturate with saline solution or sterile water and pull off slowly. In the case of catheterisation, it will be sufficient to cut the stitches with which the device was previously closed.
DEVICE REPLACEMENT AND CONTINUED TREATMENT
 The device should be replaced for the first time after 48 hours and then at least 2-3 times per week. The application of the device may be repeated for up to 30 days. However, the application frequency should be gradually reduced as the healing process is completed.
WARNINGS
 Do not expose the device to intense heat or open flames. Do not use after the expiry date. Do not reuse the device: the transmission of pathogens cannot be ruled out in the event of reuse. Do not use the device if the package is opened or

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Medical device for the treatment of penis lesions

COMPOSITION:
 Ingredients: oxygen-enriched olive oil (98%)
 Excipients: natural pear fragrance (2%)
 Substrate: polyurethane and polyester
 The device may contain traces (<0.4%) of the following allergens: CITRAL, GERANIOL, CITRONELLLOL, LINALOOL, LIMONENE.

PRODUCT DESCRIPTION
 NovoX® Andro is a medical device for the treatment of surgical, traumatic or ulcerative lesions of the penis. By creating a micro-environment which favours the activation of microcirculation but which is at the same time unfavourable to pathogen proliferation, and thanks also to the protective, barrier and soothing action of oxygen-enriched olive oil, NovoX® Andro aids the healing and correct cicatrization of lesions. NovoX® Andro is a medical device in the form of a single-use tubular dressing with a closed tip, impregnated on the inside with an oily gel. The shape of NovoX® Andro is specially designed to fit the anatomical form of the penis. NovoX® Andro is available in different sizes to suit different penis sizes. For size selection, please refer to the indications on the package.

PACKAGE CONTENTS
 1-piece box: 1 individually packaged impregnated substrate; 1 instruction manual (this leaflet)
 5-piece box: 5 individually packaged impregnated substrates; 1 instruction manual (this leaflet)
 6-piece box: 6 individually packaged impregnated substrates; 1 instruction manual (this leaflet)
 10-piece box: 10 individually packaged impregnated substrates; 1 instruction manual (this leaflet)

INTENDED USES
 NovoX® Andro is intended for the treatment of:
 • Vascular and pressure ulcers
 • Ulcers due to neoplastic diseases
 • Post-radiotherapy ulcers and lesions
 • Non-healing surgical wounds
 • Superficial, medium and deep burns
 • Surgical wounds
 • Traumatic wounds (e.g. abrasions, cuts and lacerations)
 • Superficial lesions
 NovoX® Andro is a non-sterile medical device. However, any previous contamination of the lesion is not a contraindication to the use of the device as tests have shown the intrinsic ability of the oxygen-enriched oil to create an environment unfavourable to the proliferation of pathogens commonly found on skin lesions, making it inherently protected against contamination by such pathogens.

CONTRAINDICATIONS
 NovoX® ANDRO must not be used on patients with known hypersensitivity to or who have had an allergic reaction to the device or its components. NovoX® Andro must not be used in case of previous serious skin reactions to wound dressings. In case of known hypersensitivity to dressings, NovoX® Andro must only be used further to medical assessment. NovoX® Andro should not be used when complete coverage of the lesion is not possible.

INTENDED USERS
 The device may be used by non-professional users. In the case of extensive lesions or lesions involving the deep layers of the skin, the health care personnel should provide prior explanation on how to apply the dressing correctly. Follow the instructions of the health care personnel regarding periodic check-ups.
SECONDARY DRESSINGS
 A secondary dressing may be applied in order to maintain the penis in a fixed position, if the doctor so requests. The use of elastic netting or equivalent is recommended.

NovoX® Andro may be used underneath underwear provided they do not apply too much pressure or cause excessive deformation of the device.

DEVICE APPLICATION
 Choose the NovoX® Andro size suitable for the penis size, using the size selection guide on the package. Make sure that the device covers the entire lesion area. For surgical wounds, NovoX® Andro should be applied immediately after suturing the lesion.
 Users should wash their hands thoroughly before and after use. Moreover, wear disposable gloves and use tools that are preferably sterile or otherwise subject to careful disinfection when applying the product.
 • Before applying the device, clean the lesion area as directed by health care personnel using physiological saline, sterile water or other appropriate solution. Pat dry with sterile gauze to remove excess water. Wipe off excess disinfectant if used.
 • Open the package containing the device by pulling the flaps at the markings on the package, immediately before using the device.
 • Take out the device, taking care not to touch the impregnated inside that will be in contact with the wounds.
 • The device can be cut to reduce its length if necessary, first thoroughly disinfecting the equipment used (or using sterile equipment) and taking care not to damage the substrate or the textile coverings on the inside and outside. The device's closed tip can also be cut off, if necessary in the judgement of the doctor. It is recommended to follow the doctor's instructions carefully before cutting the device. It is recommended not to cut the device unless explicitly instructed by your doctor.
 • Fit the device, taking care not to deform or crumple it excessively and ensuring that it sticks firmly to the area to be treated.
 • Wear underwear, taking care to use a panty liner to prevent any small spills from soiling the underwear.
 • As an alternative to the above, if catheterisation is present that makes it impossible to insert the device, the following steps can be taken after thorough disinfection of the tools used to cut the device (or using sterile tools):
 - Cut off the distal end to allow passage of the catheter.
 - Open the device longitudinally by cutting on one of the welds.
 - Apply the device to the penis by wrapping the device over its entire length.
 - Restore the tubular shape of the device with at least 3 or more stitches (depending on the length) so that the device is well adhered to the penis.
 - For the execution of the suture referred to in the previous point, it is recommended to:
 - Use detached stitches (minimum 3)
 - Use slowly absorbable (e.g. polyglyconate) or non-absorbable thread of the same size used for tissue suturing.

DEVICE REMOVAL
 • To remove the device, simply pull it off. If removal is difficult, saturate with saline solution or sterile water and pull off slowly. In the case of catheterisation, it will be sufficient to cut the stitches with which the device was previously closed.
DEVICE REPLACEMENT AND CONTINUED TREATMENT
 The device should be replaced for the first time after 48 hours and then at least 2-3 times per week. The application of the device may be repeated for up to 30 days. However, the application frequency should be gradually reduced as the healing process is completed.
WARNINGS
 Do not expose the device to intense heat or open flames. Do not use after the expiry date. Do not reuse the device: the transmission of pathogens cannot be ruled out in the event of reuse. Do not use the device if the package is opened or

damaged.
 Do not use the device if storage is not in accordance with these instructions.
 Do not use in conjunction with creams or medicinal products for topical use.
 Do not use if other dressings have already been applied: this could reduce the device's efficacy and/or cause local hypersensitivity reactions.
 Use on adhesive sutures may reduce the device's efficacy.
 Use on anatomical locations other than the intended ones or on parts significantly different from the physiological anatomy, or excessive deformation or crumpling of the device, may reduce the efficacy of the device and/or cause local hypersensitivity reactions.
 In the case of restoring the tubular shape by suturing, it is recommended to:
 - Take care not to deform or bend the device excessively
 - Take care not to excessively compress the penis
 - Take care that the device is well adhered to the entire lesion area
 - Care must be taken when applying the stitches, which must not affect the penis.
 A burning sensation or slight irritation immediately after application is normal and generally disappears within a short time. If the burning sensation is particularly intense or persistent, or if intense irritation appears, discontinue use and inform your doctor or paramedical staff.
 The non-impregnated outside of the device is fairly water-repellent, but the device should not be wet or immersed in water.
 Small leaks of oxygen-enriched oil from the edges around the dressing opening and from the edges of the areas where the medical device has been cut are normal.
 Oxygen-enriched oil may stain fabrics.
 The pungent smell is an intrinsic characteristic, indicating the release of reactive oxygen species. Keep out of the reach of children.
 If the product is used on minors, adult supervision is recommended.
 Promptly inform the manufacturer MOSS S.p.A. and the National Competent Authority if a serious accident occurs while using the device. Also inform the manufacturer MOSS S.p.A. of any side effects not mentioned in this leaflet.
 If you require additional information, please contact the manufacturer, MOSS S.p.A.

DISPOSAL
 Unused devices can be disposed of as unsorted waste. The device and any secondary dressings after use should be disposed of as potentially infected biological materials.
STORAGE
 Keep at a temperature below 25°C, in a cool, dry area far from sources of heat.
 In case of high room temperatures the product should be kept in the refrigerator (2-8°C).

CE 0477
 Date of last revision of product leaflet:
 Rev: MDR-0922 of 13/09/2022

	Keep at temperatures below 25°C		Product number		Lot number		Expiry date - yyyy/mm/dd
	Medical device		Unique Device Identifier		Read the instructions for use prior to use		Do not reuse
	Read warnings prior to use		Do not use if the package is damaged		Keep away from light		Keep dry

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